

ReliOn™ PRIME

Blood Glucose Monitoring System



83-05674A

Quick Reference Guide

IMPORTANT INFORMATION


- Test results below 70 mg/dL indicate low blood glucose (hypoglycemia)
- Test results above 240 mg/dL indicate high blood glucose (hyperglycemia)
- If results are below 70 mg/dL or above 240 mg/dL, repeat the test
- If results continue to measure below 70 mg/dL or above 240 mg/dL, consult your diabetes healthcare professional immediately
- If symptoms are not consistent with test results AND all the instructions have been followed carefully, contact your diabetes healthcare professional.
- Before testing, it is important to wash hands with soap and warm water. Dry hands thoroughly. If you use alcohol wipes to cleanse fingers, make sure the area is dry before lancing finger. Unwashed hands or residual alcohol may lead to inaccurate readings.
- Before testing, always make sure the meter is at operating temperature between 50–104°F (10–40°C). Treat test results obtained outside the optimum temperature range with caution.
- If the meter has been exposed to extreme temperatures, move the meter to an area within the operating range of 50–104°F (10–40°C).

PERFORMING A BLOOD GLUCOSE TEST


Step 1

Remove ReliOn™ Prime test strip from the bottle. Immediately replace the bottle cap tightly.

Insert test strip as shown in diagram. The meter will turn on automatically.



Make sure all parts of the display screen appear.



IMPORTANT!
Insert into meter this way

Step 2

Obtain a blood sample. Apply blood to the end of the test strip. Capillary action will pull the blood to the reaction site. The meter will start to count down when it detects that blood has been applied.


Do not add blood on top of the test strip, as you may receive an inaccurate result.

Discard strip if meter does not start to count down after blood drop is applied.

IMPORTANT: The minimum sample volume is 0.5 µL. Smaller samples may not start the meter or can start the meter count down and provide an "E13" error message. Ensure blood completely fills reaction site.



DO NOT add blood to the top of the strip. You may receive an abnormal (high) result.



Step 3

The test result and date/time is displayed in 7 seconds. The result will be automatically stored in memory.

Results will be shown in units of mg/dL.

Always record your results in your self-testing logbook along with other information such as insulin dosage, diet and exercise.



Step 4

When the test is done, remove the test strip. Dispose of the test strip carefully. The meter shuts off automatically when you remove the test strip.

A lancet should only be used once. After use, put the protective cap back on and dispose of in a safe manner to avoid accidental lancet stick injuries.




TROUBLESHOOTING

Refer to the ReliOn™ Prime Blood Glucose Monitoring System user instruction manual for further details.



Your blood glucose level is less than 20 mg/dL.

Repeat the test using a new test strip. If this message appears again, contact your healthcare professional immediately.



Your blood glucose level is more than 600 mg/dL.

Repeat the test using a new test strip. If this message appears again, contact your healthcare professional immediately.



The battery has been replaced.

Check that the time and date are set correctly.



A used strip has been inserted or the meter is having difficulty reading the strip.

Repeat the test with a new test strip.



The battery power is low.

There is only enough power for about 20 more tests. Replace with a new CR2032 battery as soon as possible.




The battery power is too low for testing.

Replace with a new CR2032 battery.




Thermometer flashing

The meter is outside the optimum temperature range of 50–104°F (10–40°C). Treat test result obtained where the thermometer is flashing with caution. The flashing thermometer will be displayed on all screens during glucose testing until the meter has reached acceptable temperature. Readings obtained with the flashing thermometer will show a thermometer when results are recalled in memory.



Thermometer symbol - E6.

Temperature outside of acceptable range of 41–113°F (5–45°C). You will not obtain a reading. Allow the meter to warm up or cool down slowly to between 50–104°F (10–40°C) or at least to 41–113°F (5–45°C).



Not enough blood was added to test strip.

Retest with new test strip. Apply at least 0.5 µL of blood to sample application tip.



Retest with new test strip. Apply at least 0.5 µL of blood to sample application tip.

For all other errors, call Customer Service at 855.776.0662.


PERFORMING A BLOOD GLUCOSE TEST FROM AN ALTERNATE SITE

Consider AST when:

- Testing before a meal.
- You are in a fasting state.
- Two hours have passed since a meal.
- Two hours have passed since insulin dosing.
- Two hours have passed since physical activity.


Step 1

Choose a Site: Select a soft, fleshy area on the palm that is free of visible veins and moles and away from bone. Wash the area with warm, soapy water. Rinse and dry completely. If you use alcohol wipes to cleanse the site, make sure the area is dry before lancing.



Step 2


Insert lancet and cock lancing device.




Step 3

Remove ReliOn™ Prime test strip from the bottle. Immediately replace the bottle cap tightly.

Insert test strip as shown in diagram. Be sure to align test strip with test strip port.




Make sure all parts of the display screen appear.




IMPORTANT!
Insert into meter this way

Step 4


a) **Vigorously rub the selected area 5 to 10 seconds** until it starts to feel warm to the touch.



b) Firmly hold the cocked lancing device against the clean skin for 5 to 10 seconds.



c) Press the release button on the lancing device to lance the skin. **Continue to hold the lancing device firmly against the skin** until a 0.5 µL blood drop forms.




d) Once a large enough drop of blood has formed, remove the lancing device.

IMPORTANT


Repeat blood draw if fluid is not a deep red color.

Step 5

Pick up the meter. Touch the end of the test strip (at a 90° angle) to the blood sample until the reaction site is full. Immediately remove the meter and test strip from the blood drop.



CORRECT



INCORRECT

IMPORTANT

- Was the sample deep red in color and free of watery fluid?
- Was the drop large enough (at least 0.5 µL)?
- Did the reaction site fill with enough blood all at once?

IF NOT, REPEAT TEST.



Enough Blood

Not Enough Blood

Step 6

The test result and date/time is displayed in 7 seconds. The result will be automatically stored in memory.

Results will be shown in units of mg/dL.

CONSIDER YOUR RESULT.

REPEATING ALTERNATE SITE TEST:

- IF the blood sample appeared to be diluted with clear fluid.
- IF you did not vigorously rub the test site.
- IF the test was accidentally marked as a "Control" result.
- IF your result was not consistent with how you feel.
- IF the blood drop was not large enough to fill the reaction site.



Enough Blood

Not Enough Blood

Too Much Blood

Any of the above situations can lead to an inaccurate test result.

If the repeated alternate site result is still not consistent with how you feel, confirm your blood glucose level with fingertip testing.

Step 7

Always record your results in your self-testing logbook along with other information such as insulin dosage, diet, and exercise. If you do not have a logbook, call 855.776.0662, or visit ReliOn.com to obtain a free copy.

Step 8

When the test is done, remove the test strip. Dispose of the test strip carefully. The meter shuts off automatically when you remove the test strip.

A lancet should only be used once. After use, put the protective cap back on and dispose of in a safe manner to avoid accidental lancet stick injuries.

Distributed by:
Wal-Mart Stores, Inc.
Bentonville, AR 72716
855.776.0662
ReliOn.com

